

APPROVED DISINFECTANTS

Oldelft's TEE probes should be disinfected following the instructions for use which can be downloaded from Oldelft's website: www.oldelft.com/products/downloads.

Material compatibility

Oldelft has tested the disinfectants listed in this document on chemical compatibility with the materials used in the patients applied part of our TEE probes. The patient applied part is the distal end of the probe up to the latest marker on the gastroscope shaft. Refer to this table to find out which disinfectants are compatible with each probe type. In table 1, the disinfection agents are listed that are compatible with our TEE probes.

Efficacy

Oldelft has proven that its TEE probes can be disinfected effectively using a broad range of agents. In table 1, the disinfection agents are listed that are efficacious with our TEE probes.

The efficacy claim is only valid when the instructions for use are followed.

Applicability

This document is applicable to all Oldelft's Pediatric (MiniMulti), Neonatal (MicroMulti) Adult (MultiPlane) TEE probes with Item numbers in the 19x range. For older products, please contact your sales representative.

Table 1. List of approved disinfectants

Disinfectant	Manufacturer	Active ingredient	Concentration	Compatible	Efficacious
Cidex	ASP J&J	Glutaraldehyde	2,55%	●	●
Cidex OPA	ASP J&J	Ortho-phthalaldehyde	0,55%	●	●
Gigasept FF (new)	Schulke	Dimethoxytetrahydrofuran Succindialdehyde Ethanol Methanol Alkylpolyethylenglykolpolypropylene glycol-ether, 2-ethanol, 3,6-dioxa-1-dodecanol, 3,6-dioxadodecan-1-ol, DEGHE, diethylene glycol monohexyl ether, hexyl carbitol	3,2% 11,9% 5-15% 5-10% 1-5%	●	●
Gigasept PAA concentrate	Schulke	Peracetic Acid Hydrogen peroxide Acetic Acid	5% 10-20% 10-20%	●	●
Metricide	Metrex	Glutaraldehyde	2,60%	●	●
Metricide 28	Metrex	Glutaraldehyde	2,50%	●	●
Metricide Plus 30	Metrex	Glutaraldehyde	3,40%	●	●


Rev.: A	Date: 20210819	Description: Efficacy added	Approval: RJS
Designed by: K. de Vries		Category: Approved Chemicals	Sheet: 1 of 2
		Product: List of approved disinfectants	Doc. no. 250C002 02A

APPROVED DISINFECTANTS

Disinfectant	Manufacturer	Active ingredient	Concentration	Compatible	Efficacious
Metricide OPA Plus	Metrex	Ortho-phthalaldehyde	0,60%	●	●
Perasafe	Dupont	Sodium Perborate	40-60%	●	●
Rapicide PA HLD	Medivators	Peracetic Acid Hydrogen Peroxide Trisodium Phosphate Surfactant	5% 22% 4,3% 4%	● ¹	●
Revital-Ox Resert XL HLD	Steris	Hydrogen Peroxide 2-Fluoroic Acid	1.4-2.3 ≤2.50	● ²	●
Steranios 2%, 2% N.G., 2% E.C.S.	Laboratoires Anios	Glutaraldehyde	2,00%	●	●
TD5	CS Medical	Glutaraldehyde	2,65%	●	●
Tristel Trio Wipe System	Tristel	Propan-2-OL Polymeric Biguanide Hydrochloride 5-Chloro-2-Isothiazol-3-one 2-Methyl-2H-Isothiazol-3-O Chlorinedioxide Sodiumchlorite 100%	1-10% < 1% < 1% < 1% < 1%	●	●
Wavicide 01	Medical Chemical Corp.	Glutaraldehyde	2,65%	●	●

¹ For MicroMulti and MiniMulti TEE probes. Others are in progress.

² Repeated exposure to Revital-Ox Resert XL HLD might on the long term lead to discoloration of the shaft. The product will remain electrically safe.

Rev.: A	Date: 20210819	Description: Efficacy added	Approval: RJS
Designed by: K. de Vries		Category: Approved Chemicals	Sheet: 2 of 2
		Product: List of approved disinfectants	Doc. no. 250C002 02A